



Product description

Order number 95.1007.101
Product description Tourniquet, latex-free

Product characteristics

Application venous blood collection

Purity & certification

Product category no medical device | IVD
Latex-free yes
Batched yes

Packaging

Minimum order qty. 250
Type of inner packaging case
Piece(s) / inner box 250
Piece(s) / case 2000
Piece(s) / pallet 42000
Case volume 0 cbm
Weight of case 7.01 kg
EAN of inner box 4038917263214

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN case

4038917398220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT