



Product description

Order number	14.1207
Product description	Blood Culture Adapter LongNeck, For narrow bottle neck, membrane screw cap assembled, connection: blood culture bottle, for use with (Safety) Multifly® needle, 1 piece(s)/blister, 50 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

Product characteristics

Design	membrane screw cap assembled
Suitable for	use with (Safety) Multifly® needle
Connection	blood culture bottle

Size

Diameter of opening	19 mm
Length of product	74 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	9000
Depth of box	350 mm
Width of box	110 mm
Height of box	170 mm
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0041 kg
Weight of case	2.14 kg
EAN of inner box	4038917364188
EAN case	4038917364171

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

