



### Product description

Order number	14.1209
Product description	Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled, connection: blood culture bottle, for use with (Safety) Multifly® needle, 1 piece(s)/blister, 20 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

### Product characteristics

Design	membrane adapter assembled
Suitable for	use with (Safety) Multifly® needle
Connection	blood culture bottle

### Size

Diameter of opening	38 mm
Length of product	84 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile, non-pyrogenic/endotoxin-free  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 100  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 20  
Piece (s) / outer case 100  
Piece(s) / pallet 3600  
Depth of case 578 mm  
Width of case 368 mm  
Height of case 195 mm  
Case volume 0.0415 cbm  
Weight of product 0.0079 kg  
Weight of case 1.94 kg  
EAN of inner box 4038917354776  
EAN case 4038917354769

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-10-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**