



Product description

Order number	14.1219
Product description	POCT DISPENSER, for needle-free and hygienic blood dispensing, connection: S-Monovette®, for S-Monovette®, 1 piece(s)/blister, 100 piece(s)/case, non-sterile

Product characteristics

Suitable for	S-Monovette®
Connection	S-Monovette®

Size

Length of product	49 mm
-------------------	-------

Material & colours

Colour of product	translucent
-------------------	-------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**POCT DISPENSER, for needle-free and hygienic
blood dispensing**

Page 2

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	non-sterile
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of box	268 mm
Width of box	81 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	287 mm
Height of case	172 mm
Case volume	0.0184 cbm
Weight of product	0.0015 kg
Weight of case	2.51 kg
EAN of inner box	4038917432177
EAN case	4038917432160

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP
support and is valid without signature.



SARSTEDT