



### **Product description**

Order number	14.1219
Product description	POCT DISPENSER, for needle-free and hygienic blood dispensing, connection: S-Monovette® , for S-Monovette® , 1 piece(s)/blister, 100 piece(s)/case, non-sterile

### **Product characteristics**

Suitable for	S-Monovette®
Connection	S-Monovette®

### **Size**

Length of product	49 mm
-------------------	-------

### **Material & colours**

Colour of product	translucent
-------------------	-------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category	no medical device   IVD
Purity standard	non-sterile
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of inner packaging	blister
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	80000
Depth of box	268 mm
Width of box	81 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	178 mm
Height of case	287 mm
Case volume	0.0191 cbm
Weight of product	0.0015 kg
Weight of case	2.51 kg
EAN of inner box	4038917432177
EAN case	4038917432160

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

