



Product Specification

**Safety-Multifly® BKF set 21G, 21G x 3/4", green,
tube length: 200 mm, 1 piece(s)/blister**

Page 1



Product description

Order number

51.1638.945

Product description

Safety-Multifly® BKF set 21G, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 200 mm, with blood culture adapter set, set components: 1 Safety-Multifly® needle 21G with long tubing, DEHP-free (REF 85.1638.235) and 1 Blood Culture Adapter Universal (REF 14.1209), 1 piece(s)/blister, sterile, 1 piece(s)/case

Product characteristics

Design

blood culture adapter set

Coding colour

green

Safety product

yes

Connection

blood culture bottle

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Safety-Multifly® BKF set 21G, 21G x 3/4", green,
tube length: 200 mm, 1 piece(s)/blister

Page 2

Size

Diameter of needle (external)	21 G / 0,8 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	200 mm

Purity & certification

Product category	Medical device (MD)
CE certified	no CE mark
Purity standard	sterile
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	2100
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.0118 kg
Weight of case	1.42 kg
EAN of inner box	4038917357470
EAN case	4038917357463

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT