




## Product Specification

**S-Monovette® safety needle, 20G x 1 1/2", yellow,  
1 piece(s)/blister**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 85.1160.200  |
| Product description | S-Monovette® safety needle, needle (ØxL): 20G x 1 1/2", coding colour: yellow, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case |
| PZN number          | 6796224  |

### Product characteristics

|                |                    |
|----------------|--------------------|
| Coding colour  | yellow             |
| Safety product | yes                |
| Connection     | membrane screw cap |

### Size

|                               |               |
|-------------------------------|---------------|
| Diameter of needle (external) | 20 G / 0,9 mm |
| Length of needle              | 38 mm         |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® safety needle, 20G x 1 1/2", yellow,  
1 piece(s)/blister**

**Page 2**

### Purity & certification

|                  |                                       |
|------------------|---------------------------------------|
| Product category | Medical device                        |
| CE certified     | CE 0124 (DEKRA)                       |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation    | Gamma irradiation                     |
| Pyrogen-free     | yes                                   |
| Batched          | yes                                   |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 39000         |
| Depth of box                  | 268 mm        |
| Width of box                  | 111 mm        |
| Height of box                 | 71 mm         |
| Depth of case                 | 373 mm        |
| Width of case                 | 233 mm        |
| Height of case                | 290 mm        |
| Case volume                   | 0.0252 cbm    |
| Weight of product             | 0.0026 kg     |
| Weight of case                | 2.48 kg       |
| EAN of inner box              | 4038917126021 |
| EAN case                      | 4038917361910 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**