



Product description

Order number	85.1372.200
Product description	S-Monovette® safety needle, needle (ØxL): 20G x 1", coding colour: yellow, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free

Product characteristics

Coding colour	yellow
Safety product	yes
Connection	membrane screw cap

Size

Diameter of needle (external)	20 G / 0,9 mm
Length of needle	25,4 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**S-Monovette® safety needle, 20G x 1", yellow, 1
piece(s)/blister**

Page 2

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0124 (Dekra)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Lot or serial number	yes

Packaging

Minimum order qty.	500
Type of inner packaging	blister
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	39000
Depth of box	268 mm
Width of box	111 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	233 mm
Height of case	290 mm
Case volume	0,0252 cbm
Weight of product	0,0026 kg
Weight of case	2,47 kg
EAN of inner box	4038917127943
EAN case	4038917361941

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-06-23

This document was prepared by EDP support and is valid without signature.



SARSTEDT