




Product Specification
**S-Monovette® needle, 21G x 1", green, 1
piece(s)/blister**

Page 1



Product description

Order number	85.1373
Product description	S-Monovette® needle, needle (ØxL): 21G x 1", coding colour: green, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case

Product characteristics

Coding colour	green
Safety product	no
Connection	membrane screw cap

Size

Diameter of needle (external)	21 G / 0,8 mm
Length of needle	25.4 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
S-Monovette® needle, 21G x 1", green, 1
piece(s)/blister

Page 2

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	78000
Depth of box	268 mm
Width of box	111 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	233 mm
Height of case	290 mm
Case volume	0.0252 cbm
Weight of product	0.0017 kg
Weight of case	2.92 kg
EAN of inner box	4038917086387
EAN case	4038917004930

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT