



Product description

Order number	85.1440
Product description	S-Monovette® needle, needle (ØxL): 22G x 1 1/2", coding colour: black, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case

Product characteristics

Coding colour	black
Safety product	no
Connection	membrane screw cap

Size

Diameter of needle (external)	22 G / 0,7 mm
Length of needle	38 mm

Purity & certification

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® needle, 22G x 1 1/2", black

Page 2

Batched yes

Packaging

Minimum order qty.	1000
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	78000
Depth of box	268 mm
Width of box	111 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	233 mm
Height of case	290 mm
Case volume	0.0252 cbm
Weight of product	0.0017 kg
Weight of case	2.92 kg
EAN of inner box	4038917086653
EAN case	4038917005180

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-01-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT