



Product description

Order number	85.1440
Product description	S-Monovette® needle, needle (ØxL): 22G x 1 1/2", coding colour: black, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case

Product characteristics

Coding colour	black
Safety product	no
Connection	membrane screw cap

Size

Diameter of needle (external)	22 G / 0,7 mm
Length of needle	38 mm

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



