



Product description

Order number	85.1638.035
Product description	Multifly® needle, needle (ØxL): 21 G x 3/4", coding colour: green, tube length: 200 mm, with Multi-adapter, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case

Product characteristics

Design	Multi-adapter
Coding colour	green
Safety product	no
Connection	membrane screw cap

Size

Diameter of needle (external)	21 G / 0,8 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	200 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Multifly® needle, 21G x 3/4", green, tube length:
200 mm, 1 piece(s)/blister

Page 2

Purity & certification

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

Packaging

Minimum order qty.	480
Type of smallest subpackaging	blister
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	20640
Depth of box	224 mm
Width of box	214 mm
Height of box	171 mm
Depth of case	442 mm
Width of case	358 mm
Height of case	240 mm
Case volume	0.038 cbm
Weight of product	0.0032 kg
Weight of case	3.3 kg
EAN of inner box	4038917423779
EAN case	4038917423489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-01-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT