




Product Specification
Safety-Multifly® needle, 23G x 3/4", blue, tube
length: 200 mm, 1 piece(s)/blister

Page 1



Product description

Order number	85.1640.200
Product description	Safety-Multifly® needle, needle (ØxL): 23G x 3/4", coding colour: blue, tube length: 200 mm, with multi-adapter, screw-on guide sleeve, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case
PZN number	6797436

Product characteristics

Design	multi-adapter, screw-on guide sleeve
Coding colour	blue
Safety product	yes
Connection	membrane screw cap

Size

Diameter of needle (external)	23 G / 0,6 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	200 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Safety-Multifly® needle, 23G x 3/4", blue, tube
length: 200 mm, 1 piece(s)/blister

Page 2

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

Packaging

Minimum order qty.	480
Type of smallest subpackaging	blister
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	20640
Depth of box	224 mm
Width of box	214 mm
Height of box	171 mm
Depth of case	442 mm
Width of case	358 mm
Height of case	240 mm
Case volume	0.038 cbm
Weight of product	0.004 kg
Weight of case	3.64 kg
EAN of inner box	4038917423939
EAN case	4038917423649

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT