



Product Specification
Safety-Multifly® needle, 23G x 3/4", blue, tube
length: 240 mm, 1 piece(s)/blister



Product description

Order number	85.1640.201
Product description	Safety-Multifly® needle, needle (ØxL): 23G x 3/4", coding colour: blue, tube length: 240 mm, with Luer (f), DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case

Product characteristics

Design	Luer (f)
Coding colour	blue
Safety product	yes
Connection	Luer (f)

Size

Diameter of needle (external)	23 G / 0,6 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	240 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Safety-Multifly® needle, 23G x 3/4", blue, tube
length: 240 mm, 1 piece(s)/blister

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

Packaging

Minimum order qty.	480
Type of inner packaging	blister
Piece(s) / inner box	120
Piece(s) / case	480
Piece(s) / pallet	20640
Depth of box	224 mm
Width of box	214 mm
Height of box	171 mm
Depth of case	442 mm
Width of case	358 mm
Height of case	240 mm
Case volume	0.038 cbm
Weight of product	0.003 kg
Weight of case	3.32 kg
EAN of inner box	4038917423946
EAN case	4038917423656

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

