




**Product Specification**  
**Safety-Multifly® needle, 23G x 3/4", blue, tube**  
**length: 80 mm, 1 piece(s)/blister**

**Page 1**



**Product description**

Order number	85.1640.205
Product description	Safety-Multifly® needle, needle (ØxL): 23G x 3/4", coding colour: blue, tube length: 80 mm, with Multi-adapter, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case
PZN number	6797532

**Product characteristics**

Design	Multi-adapter
Coding colour	blue
Safety product	yes
Connection	membrane screw cap

**Size**

Diameter of needle (external)	23 G / 0,6 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	80 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Safety-Multifly® needle, 23G x 3/4", blue, tube**  
**length: 80 mm, 1 piece(s)/blister**

**Page 2**

**Purity & certification**

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

**Packaging**

Minimum order qty.	480
Type of smallest subpackaging	blister
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	20640
Depth of box	224 mm
Width of box	214 mm
Height of box	171 mm
Depth of case	442 mm
Width of case	358 mm
Height of case	240 mm
Case volume	0.038 cbm
Weight of product	0.0033 kg
Weight of case	3.31 kg
EAN of inner box	4038917423960
EAN case	4038917423670

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**