



### Product description

Order number	85.925
Product description	Micro needle, for collecting very small amounts of blood, needle (ØxL): 25G x 3/4", coding colour: orange, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case

### Product characteristics

Application	for collecting very small amounts of blood
Coding colour	orange
Safety product	no

### Size

Diameter of needle (external)	25 G / 0,5 mm
Length of needle	19 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

### **Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	55000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	230 mm
Height of case	160 mm
Case volume	0.0124 cbm
Weight of product	0.0013 kg
Weight of case	1.36 kg
EAN of inner box	4038917281843
EAN case	4038917281836

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

