



Product description

Order number	26.323
Product description	Sample tube, nominal volume: 10 ml, cap white, preparation: Clotting activator, push cap, (LxØ): 95 x 16.8 mm, round base, 100 piece(s)/bag

Product characteristics

Type of preparation	Clotting activator
Colour code	EU
Graduation	no
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, Serum CAT, 10 ml, cap white, (LxØ):
95 x 16.8 mm

Size

Sample volume	10 ml
Diameter	16.8 mm
Length of product	95 mm
Length excluding cap	95 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white
Tube material	Polypropylene (PP)

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0066 kg
Weight of case	7.27 kg
EAN of inner box	4038917081146
EAN case	4038917000284

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

