



## Product Specification

Sample tube, Lithium heparin LH, 10 ml, cap orange, (LxØ): 95 x 16.8 mm, with print

Page 1



### Product description

Order number	26.369
Product description	Sample tube, nominal volume: 10 ml, cap orange, preparation: Lithium heparin LH, push cap, (LxØ): 95 x 16.8 mm, round base, with print, 100 piece(s)/bag

### Product characteristics

Type of preparation	Lithium heparin LH
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with print
Colour of print/label	orange
Graduation	yes
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Sample tube, Lithium heparin LH, 10 ml, cap orange, (LxØ): 95 x 16.8 mm, with print**

**Page 2**

### Size

Sample volume	10 ml
Diameter	16.8 mm
Length of product	95 mm
Length excluding cap	95 mm

### Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	orange
Tube material	Polypropylene (PP)

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0066 kg
Weight of case	7.28 kg
EAN of inner box	4038917081344
EAN case	4038917000482

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**