



**Product Specification**  
**Sample tube, Serum CAT, 10 ml, cap white, (LxØ):**  
**100 x 15.7 mm**

**Page 1**



**Product description**

Order number	28.364
Product description	Sample tube, nominal volume: 10 ml, cap white, preparation: Clotting activator, push cap, (LxØ): 100 x 15.7 mm, round base, 500 piece(s)/bag

**Product characteristics**

Type of preparation	Clotting activator
Colour code	EU
Graduation	no
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, Serum CAT, 10 ml, cap white, (LxØ):**  
**100 x 15.7 mm**

**Size**

Sample volume	10 ml
Diameter	15.7 mm
Length of product	100 mm
Length excluding cap	100 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	310 mm
Case volume	0.0491 cbm
Weight of product	0.0069 kg
Weight of case	7.58 kg
EAN of inner box	4038917081320
EAN case	4038917000468

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

