



## Product description

Order number	29.373.004
Product description	Sample tube, nominal volume: 7 ml, cap black, preparation: Clotting activator, push cap, (LxØ): 100 x 13 mm, round base, with print, 100 piece(s)/bag

## Product characteristics

Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with print
Colour of print/label	white
Graduation	no
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, Serum CAT, 7 ml, cap black, (LxØ):**  
**100 x 13 mm, with print**

**Page 2**

**Size**

Sample volume	7 ml
Diameter	13 mm
Length of product	100 mm
Length excluding cap	100 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	black
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of case	375 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0354 cbm
Weight of product	0.0047 kg
Weight of case	5.26 kg
EAN of inner box	4038917105859
EAN case	4038917335973

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

