



**Product Specification**  
**Sample tube, Lithium heparin LH, 4.5 ml, cap green, (LxØ): 75 x 13 mm, with paper label**

**Page 1**



**Product description**

Order number	32.331.100
Product description	Sample tube, nominal volume: 4.5 ml, cap green, preparation: Lithium heparin, push cap, (LxØ): 75 x 13 mm, round base, with paper label, 100 piece(s)/bag

**Product characteristics**

Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Sample tube, Lithium heparin LH, 4.5 ml, cap green, (LxØ): 75 x 13 mm, with paper label**

**Page 2**

### Size

Sample volume	4.5 ml
Diameter	13 mm
Length of product	75 mm
Length excluding cap	75 mm

### Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0039 kg
Weight of case	4.37 kg
EAN of inner box	4038917105439
EAN case	4038917010573

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**