



Product Specification

**Sample tube, Lithium heparin LH, 4 ml, cap green,
(LxØ): 75 x 12 mm, with paper label**

Page 1



Product description

Order number	36.336.100
Product description	Sample tube, nominal volume: 4 ml, cap green, preparation: Lithium heparin LH, push cap, (LxØ): 75 x 12 mm, round base, with paper label, 100 piece(s)/bag

Product characteristics

Type of preparation	Lithium heparin LH
Preparation concentration	16 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Sample tube, Lithium heparin LH, 4 ml, cap green,
(LxØ): 75 x 12 mm, with paper label**

Page 2

Size

Sample volume	4 ml
Diameter	12 mm
Length of product	75 mm
Length excluding cap	75 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0035 kg
Weight of case	7.69 kg
EAN of inner box	4038917105521
EAN case	4038917010641

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT