



Product Specification
Sample tube, EDTA K3E, 2 ml, cap red, (LxØ): 55 x 12 mm, with print

Page 1



Product description

Order number	38.346
Product description	Sample tube, nominal volume: 2 ml, cap red, preparation: K3 EDTA, push cap, (LxØ): 55 x 12 mm, round base, with print, 100 piece(s)/bag

Product characteristics

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with print
Colour of print/label	red
Graduation	yes
Cap	Skirted push cap
Closure type	push cap
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Sample tube, EDTA K3E, 2 ml, cap red, (LxØ): 55 x 12 mm, with print

Page 2

Size

Sample volume	2 ml
Diameter	12 mm
Length of product	55 mm
Length excluding cap	55 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0026 kg
Weight of case	5.81 kg
EAN of inner box	4038917081283
EAN case	4038917000420

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

