



Product description

| | |
|---------------------|---|
| Order number | 41.1392.105 |
| Product description | Micro sample tube Serum CAT, venous blood collection, preparation: Clotting activator, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: red, colour code ISO, with paper label, label/print: red, skirted conical base, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|-------------------------|
| Type of collection | venous |
| Type of preparation | Clotting activator |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | red |
| Graduation | yes |
| Closure type | screw cap |
| Application | venous blood collection |
| Base shape | skirted conical base |
| Volume of work | 1,3 ml |
| Recommended centrifugation | 10000 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Micro sample tube Serum CAT, 1.3 ml, screw cap,
ISO

Page 2

Centrifugation time 5 Min

Size

Sample volume 1.3 ml
Diameter 10.8 mm
Length of product 47 mm
Length including cap 47 mm
Length excluding cap 45 mm

Material & colours

Colour of product transparent
Closure material Polypropylene (PP)
Colour of cap red
Tube material Polypropylene (PP)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 96000
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.0016 kg
Weight of case 1.94 kg
EAN of inner box 4038917128285
EAN case 4038917021630

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT