



## Product description

Order number	41.1395.105
Product description	Micro sample tube EDTA K3E, venous blood collection, preparation: K3 EDTA, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: violet, colour code ISO, with paper label, label/print: violet, skirted conical base, 100 piece(s)/bag

## Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	violet
Graduation	yes
Closure type	screw cap
Application	venous blood collection
Base shape	skirted conical base
Volume of work	1,3 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	1.3 ml
Diameter	10.8 mm
Length of product	47 mm
Length including cap	47 mm
Length excluding cap	45 mm

### **Material & colours**

Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	violet
Tube material	Polypropylene (PP)

### **Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0016 kg
Weight of case	1.9 kg
EAN of inner box	4038917128407
EAN case	4038917021708

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

