



## Product description

Order number	41.1500.005
Product description	Micro sample tube Serum Gel CAT, venous blood collection, preparation: Clotting activator/gel, nominal volume: 1.1 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: brown, colour code EU/ISO, with plastic label, label/print: brown, skirted conical base, 100 piece(s)/bag

## Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	brown
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base
Volume of work	1,1 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube Serum Gel CAT, 1.1 ml, push cap, EU/ISO**

**Page 2**

Recommended centrifugation 10000 U/min  
Centrifugation time 5 Min

**Size**

Sample volume 1.1 ml  
Diameter 10.8 mm  
Length of product 44 mm  
Length including cap 44 mm  
Length excluding cap 40 mm

**Material & colours**

Colour of product transparent  
Closure material Low Density Polyethylene (LD-PE)  
Colour of cap brown  
Tube material Polypropylene (PP)

**Purity & certification**

Product category In-vitro diagnostic (IVD)  
CE certified CE - manufacturer's self-declaration  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 96000  
Depth of case 348 mm  
Width of case 278 mm  
Height of case 140 mm  
Case volume 0.0135 cbm  
Weight of product 0.0018 kg  
Weight of case 2.11 kg  
EAN of inner box 4038917129473  
EAN case 4038917022248

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

