



Product description

Order number	41.1501.005
Product description	Micro sample tube Serum, venous blood collection, preparation: clotting activator, nominal volume: 1.3 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: white, colour code EU, with plastic label, label/print: white, skirted conical base, 100 piece(s)/bag

Product characteristics

Type of collection	venous
Type of preparation	clotting activator
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	white
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base
Volume of work	1,3 ml
Recommended centrifugation	10000 U/min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Centrifugation time 5 Min

Size

Sample volume 1.3 ml
Diameter 10.8 mm
Length of product 44 mm
Length including cap 44 mm
Length excluding cap 40 mm

Material & colours

Colour of product transparent
Closure material Low Density Polyethylene (LD-PE)
Colour of cap white
Tube material Polypropylene (PP)

Purity & certification

Product category In-vitro diagnostic, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 96000
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.0016 kg
Weight of case 1.87 kg
EAN of inner box 4038917129497
EAN case 4038917022255

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

