



Product Specification
**Micro sample tube Lithium heparin LH, 1.3 ml,
push cap, EU**

Page 1



Product description

Order number	41.1503.015
Product description	Micro sample tube Lithium heparin LH, venous blood collection, preparation: Lithium heparin, nominal volume: 1.3 ml, Ø: 10.8 mm, push cap, cap: orange, colour code EU, with paper label, label/print: orange, skirted conical base, 100 piece(s)/bag

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	orange
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Micro sample tube Lithium heparin LH, 1.3 ml,
push cap, EU**

Page 2

Volume of work 1,3 ml
Recommended centrifugation 2000 x g
Centrifugation time 5 Min

Size

Sample volume 1.3 ml
Diameter 10.8 mm
Length of product 44 mm
Length excluding cap 40 mm

Material & colours

Colour of product transparent
Closure material Low Density Polyethylene (LD-PE)
Colour of cap orange
Tube material Polypropylene (PP)

Purity & certification

Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 96000
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.0016 kg
Weight of case 1.89 kg
EAN of inner box 4038917310192
EAN case 4038917310185

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

