



### Product description

Order number	47.556.001
Product description	Sample tube, nominal volume: 3 ml, cap red, preparation: K3 EDTA, screw cap, (LxØ): 82 x 11.5 mm, skirted conical base, with paper label, 100 piece(s)/bag

### Product characteristics

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	red
Graduation	fill mark
Closure type	screw cap
Base shape	skirted conical base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Sample tube, EDTA K3E, 3 ml, cap red, (LxØ): 82 x 11.5 mm, with paper label**

**Page 2**

### Size

Sample volume	3 ml
Diameter	11.5 mm
Length of product	82 mm
Length excluding cap	66 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0028 kg
Weight of case	6.37 kg
EAN of inner box	4038917113274
EAN case	4038917013956

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**