



**Product Specification**  
**Sample tube, EDTA K3E, 3 ml, cap red, (LxØ): 82 x 11.5 mm, with print**

**Page 1**



**Product description**

Order number 47.556  
Product description Sample tube, nominal volume: 3 ml, cap red, preparation: K3 EDTA, screw cap, (LxØ): 82 x 11.5 mm, flat base, with print, 100 piece(s)/bag

**Product characteristics**

Type of preparation K3 EDTA  
Preparation concentration 1.6 mg/ml blood  
Colour code EU  
Label/ Print with print  
Colour of print/label red  
Graduation yes  
Closure type screw cap  
Base shape flat base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, EDTA K3E, 3 ml, cap red, (LxØ): 82 x 11.5 mm, with print**

**Page 2**

**Size**

Sample volume	3 ml
Diameter	11.5 mm
Length of product	82 mm
Length including cap	82 µm
Length excluding cap	66 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0027 kg
Weight of case	6.2 kg
EAN of inner box	4038917082440
EAN case	4038917001533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**