



Product Specification
Sample tube, EDTA K3E, 3 ml, cap red, (LxØ): 82 x 11.5 mm, with print

Page 1



Product description

| | |
|---------------------|---|
| Order number | 47.556 |
| Product description | Sample tube, nominal volume: 3 ml, cap red, preparation: K3 EDTA, screw cap, (LxØ): 82 x 11.5 mm, flat base, with print, 100 piece(s)/bag |

Product characteristics

| | |
|---------------------------|-----------------|
| Type of preparation | K3 EDTA |
| Preparation concentration | 1.6 mg/ml blood |
| Colour code | EU |
| Label/ Print | with print |
| Colour of print/label | red |
| Graduation | yes |
| Closure type | screw cap |
| Base shape | flat base |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, EDTA K3E, 3 ml, cap red, (LxØ): 82 x 11.5 mm, with print

Page 2

Size

| | |
|----------------------|---------|
| Sample volume | 3 ml |
| Diameter | 11.5 mm |
| Length of product | 82 mm |
| Length including cap | 82 µm |
| Length excluding cap | 66 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | red |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 72000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 300 mm |
| Case volume | 0.0475 cbm |
| Weight of product | 0.0027 kg |
| Weight of case | 6.2 kg |
| EAN of inner box | 4038917082440 |
| EAN case | 4038917001533 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT