



## Product Specification

Sample tube, Lithium heparin LH, 5 ml, cap orange, (LxØ): 57 x 16.5 mm, with paper label

Page 1



### Product description

Order number	49.354.001
Product description	Sample tube, nominal volume: 5 ml, cap orange, preparation: Heparin, screw cap, (LxØ): 57 x 16.5 mm, flat base, with paper label, 100 piece(s)/bag

### Product characteristics

Type of preparation	Heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	orange
Graduation	fill mark
Closure type	screw cap
Base shape	flat base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Sample tube, Lithium heparin LH, 5 ml, cap orange, (LxØ): 57 x 16.5 mm, with paper label**

**Page 2**

### Size

Sample volume	5 ml
Diameter	16.5 mm
Length of product	57 mm
Length excluding cap	57 mm

### Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	orange
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of case	375 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0354 cbm
Weight of product	0.0049 kg
Weight of case	5.43 kg
EAN of inner box	4038917105613
EAN case	4038917010719

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**