



Product Specification

S-Monovette® Metal analysis LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with paper label

Page 1



Product description

Order number	01.1604.400
Product description	S-Monovette® Metal analysis LH, preparation: Lithium-Heparin, liquid, 7.5 ml, membrane screw cap, cap orange, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium-Heparin, liquid
Preparation concentration	19 IU/ml blood
Label/ Print	with paper label
Colour of print/label	orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® Metal analysis LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with paper label

Page 2

Connection membrane screw cap

Size

Sample volume	7.5 ml
Diameter	15 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	32000
Depth of box	189 mm
Width of box	134 mm
Height of box	82 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Metal analysis LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with paper label

Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0073 kg
Weight of case	4.49 kg
EAN of inner box	4038917131926
EAN case	4038917024044

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

