



**Product Specification**  
**S-Monovette® EDTA K3E, 7.5 ml, cap violet,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 1**



**Product description**

|                     |                                                                                                                                                                                                       |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 01.1605.100                                                                                                                                                                                           |
| Product description | S-Monovette® EDTA K3E, preparation: K3 EDTA, 7.5 ml, membrane screw cap, cap violet, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: violet, 50 piece(s)/case, sterile |

**Product characteristics**

|                           |                    |
|---------------------------|--------------------|
| Type of collection        | venous             |
| Type of preparation       | K3 EDTA            |
| Preparation concentration | 1.6 mg/ml blood    |
| Colour code               | ISO                |
| Label/ Print              | with paper label   |
| Colour of print/label     | violet             |
| Graduation                | fill mark          |
| Cap                       | membrane screw cap |
| Closure type              | screw cap          |
| Liquid-tight              | yes                |
| Connection                | membrane screw cap |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® EDTA K3E, 7.5 ml, cap violet,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 2**

**Size**

|                                  |        |
|----------------------------------|--------|
| Sample volume                    | 7.5 ml |
| Outer diameter of cap            | 18 mm  |
| Diameter                         | 15 mm  |
| Length including cap             | 108 mm |
| Length excluding cap             | 92 mm  |
| Length including cap and plunger | 129 mm |

**Material & colours**

|                         |                                   |
|-------------------------|-----------------------------------|
| Colour of product       | transparent                       |
| Closure material        | High Density Polyethylene (HD-PE) |
| Colour of cap           | violet                            |
| Membrane material       | rubber                            |
| Tube material           | Polypropylene (PP)                |
| Piston material         | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS)                  |

**Purity & certification**

|                           |                                      |
|---------------------------|--------------------------------------|
| Satisfies the requirement | IATA, ADR                            |
| Product category          | In-vitro diagnostic (IVD)            |
| CE certified              | CE - manufacturer's self-declaration |
| Purity standard           | sterile                              |
| Sterilisation             | Electron irradiation                 |
| Batched                   | yes                                  |

**Packaging**

|                               |        |
|-------------------------------|--------|
| Minimum order qty.            | 500    |
| Type of smallest subpackaging | case   |
| Piece(s) / inner box          | 50     |
| Piece (s) / outer case        | 500    |
| Piece(s) / pallet             | 32000  |
| Depth of box                  | 189 mm |
| Width of box                  | 134 mm |
| Height of box                 | 82 mm  |
| Depth of case                 | 440 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® EDTA K3E, 7.5 ml, cap violet,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 3**

|                   |               |
|-------------------|---------------|
| Width of case     | 400 mm        |
| Height of case    | 144 mm        |
| Case volume       | 0.0253 cbm    |
| Weight of product | 0.0077 kg     |
| Weight of case    | 4.67 kg       |
| EAN of inner box  | 4038917132039 |
| EAN case          | 4038917024143 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**