



## Product Specification

**S-Monovette® Lithium heparin LH, fluid, 7.5 ml,  
cap green, (LxØ): 92 x 15 mm, with paper label**

**Page 1**



### Product description

Order number	01.1608.100
Product description	S-Monovette® Lithium heparin LH, fluid, preparation: Lithium-Heparin, liquid, 7.5 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium-Heparin, liquid
Preparation concentration	19 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Lithium heparin LH, fluid, 7.5 ml,  
cap green, (LxØ): 92 x 15 mm, with paper label**

**Page 2**

Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	7.5 ml
Diameter	15 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	32000
Depth of box	189 mm
Width of box	134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-21

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**



**Product Specification**

**S-Monovette® Lithium heparin LH, fluid, 7.5 ml,  
cap green, (LxØ): 92 x 15 mm, with paper label**

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0076 kg
Weight of case	4.64 kg
EAN of inner box	4038917132084
EAN case	4038917024181

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

