




Product Specification

**S-Monovette® Metal analysis LH, 7.5 ml, cap blue,
(LxØ): 92 x 15 mm, with paper label**

Page 1



Product description

Order number	01.1609.201
Product description	S-Monovette® Metal analysis LH, preparation: Lithium-Heparin, liquid, 7.5 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium-Heparin, liquid
Preparation concentration	19 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® Metal analysis LH, 7.5 ml, cap blue,
(LxØ): 92 x 15 mm, with paper label**

Width of box	134 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0087 kg
Weight of case	4.41 kg
EAN of inner box	4038917531689
EAN case	4038917531672

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

