



Product Specification
S-Monovette® EDTA Gel K2E, 7.5 ml, cap red,
(LxØ): 92 x 15 mm, with paper label

Page 1



Product description

Order number	01.1621.001
Product description	S-Monovette® EDTA Gel K2E, preparation: K2 EDTA Gel, 7.5 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: red/brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	K2 EDTA Gel
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	red/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® EDTA Gel K2E, 7.5 ml, cap red,
(LxØ): 92 x 15 mm, with paper label

Page 3

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0087 kg
Weight of case	5.21 kg
EAN of inner box	4038917312707
EAN case	4038917312691

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT