



Product Specification

**S-Monovette® neutral Z, 7.5 ml, cap white, (LxØ):
92 x 15 mm, with paper label pre-barcoded**

Page 1

Product description

Order number	01.1728.088
Product description	S-Monovette® neutral Z, without preparation, 7.5 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label pre-barcoded, label/print: black, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	without preparation
Colour code	ISO
Label/ Print	with paper label pre-barcoded
Colour of print/label	black
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	7.5 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® neutral Z, 7.5 ml, cap white, (LxØ):
92 x 15 mm, with paper label pre-barcoded**

Page 2

Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	32000
Depth of box	189 mm
Width of box	134 mm
Height of box	82 mm
Depth of case	435 mm
Width of case	400 mm
Height of case	145 mm
Case volume	0.0252 cbm
Weight of product	0.0076 kg
Weight of case	4.63 kg
EAN of inner box	4066916033453
EAN case	4066916033460

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT