



Product description

Order number	01.2040.001
Product description	S-Monovette® cfDNA Exact, preparation: cfDNA stabiliser, 9.2 ml, membrane screw cap, cap raspberry coloured, (LxØ) without cap: 100 x 15 mm, with paper label, label/print: dark red, 20 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	cfDNA stabiliser
Label/ Print	with paper label
Colour of print/label	dark red
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® cfDNA Exact, 9.2 ml, cap raspberry coloured, (LxØ): 100 x 15 mm, with paper label

Page 2

Size

Sample volume	9.2 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	116 mm
Length excluding cap	100 mm
Length including cap and plunger	136 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	raspberry coloured
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	80
Type of smallest subpackaging	case
Piece(s) / inner box	20
Piece (s) / outer case	80
Piece(s) / pallet	8640
Depth of box	184 mm
Width of box	141 mm
Height of box	83 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® cfDNA Exact, 9.2 ml, cap raspberry coloured, (LxØ): 100 x 15 mm, with paper label

Page 3

Depth of case	388 mm
Width of case	345 mm
Height of case	105 mm
Case volume	0.0141 cbm
Weight of product	0.0099 kg
Weight of case	1.9 kg
EAN of inner box	4038917576666
EAN case	4038917576659

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT