



Product Specification

S-Monovette® Lithium heparin LH, 9 ml, cap orange, (LxØ): 92 x 16 mm, with plastic label

Page 1



Product description

Order number	02.1065
Product description	S-Monovette® Lithium heparin LH, LH, preparation: Lithium heparin, 9 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 92 x 16 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Lithium heparin LH, 9 ml, cap orange, (LxØ): 92 x 16 mm, with plastic label

Page 2

Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

Size

Sample volume	9 ml
Outer diameter of cap	19 mm
Diameter	16 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	204 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® Lithium heparin LH, 9 ml, cap orange, (LxØ): 92 x 16 mm, with plastic label

Width of box	134 mm
Height of box	88 mm
Depth of case	648 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0088 kg
Weight of case	5.33 kg
EAN of inner box	4038917085236
EAN case	4038917003926

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

