



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 10 ml,
cap green, (LxØ): 92 x 16 mm, with paper label**

Page 1



Product description

Order number	02.1067.001
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 9NC, preparation: Citrate 9NC, 10 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 92 x 16 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 10 ml, cap green, (LxØ): 92 x 16 mm, with paper label

Page 2

Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

Size

Sample volume	10 ml
Diameter	16 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	204 mm
Width of box	134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 10 ml,
cap green, (LxØ): 92 x 16 mm, with paper label**

Page 3

Height of box	88 mm
Depth of case	648 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0088 kg
Weight of case	5.37 kg
EAN of inner box	4038917125239
EAN case	4038917019859

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP
support and is valid without signature.



SARSTEDT