



## Product description

Order number	02.1726.001
Product description	S-Monovette® neutral Z, without preparation, 9 ml, membrane screw cap, cap white, (LxØ) without cap: 92 x 16 mm, with paper label, label/print: black, 50 piece(s)/case, sterile

## Product characteristics

Type of collection	venous
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	black
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® neutral Z, 9 ml, cap white, (LxØ): 92 x 16 mm, with paper label**

**Page 2**

### Size

Sample volume	9 ml
Diameter	16 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	204 mm
Width of box	134 mm
Height of box	88 mm
Depth of case	648 mm
Width of case	303 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® neutral Z, 9 ml, cap white, (LxØ): 92  
x 16 mm, with paper label**

Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0078 kg
Weight of case	4.84 kg
EAN of inner box	4038917133203
EAN case	4038917024723

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

