




Product Specification
S-Monovette® Serum CAT, 5.5 ml, cap white,
(LxØ): 75 x 15 mm, with paper label

Page 1



Product description

Order number	03.1397.001
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 5.5 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: black, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	black
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 5.5 ml, cap white,
(LxØ): 75 x 15 mm, with paper label

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 5.5 ml
Diameter 15 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap white
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 34000
Depth of box 189 mm
Width of box 117 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 5.5 ml, cap white,
(LxØ): 75 x 15 mm, with paper label

Page 3

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.007 kg
Weight of case	4.28 kg
EAN of inner box	4038917128438
EAN case	4038917021722

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT