



Product Specification

S-Monovette® Lithium heparin gel LH, 4.7 ml, cap orange, (LxØ): 75 x 15 mm, with plastic label

Page 1



Product description

| | |
|---------------------|---|
| Order number | 03.1631 |
| Product description | S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 4.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 15 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|---------------------|
| Type of collection | venous |
| Type of preparation | Lithium heparin gel |
| Preparation concentration | 25 IU/ml blood |
| Colour code | EU |
| Label/ Print | with plastic label |
| Colour of print/label | transparent/brown |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2500 x g |
| Centrifugation temperature | 20 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Lithium heparin gel LH, 4.7 ml, cap orange, (LxØ): 75 x 15 mm, with plastic label

Page 2

| | |
|---------------------|--------------------|
| Centrifugation time | 15 Min |
| Liquid-tight | yes |
| Connection | membrane screw cap |

Size

| | |
|----------------------------------|--------|
| Sample volume | 4.7 ml |
| Diameter | 15 mm |
| Length including cap | 91 mm |
| Length excluding cap | 75 mm |
| Length including cap and plunger | 112 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | orange |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 26000 |
| Depth of box | 189 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin gel LH, 4.7 ml, cap orange, (LxØ): 75 x 15 mm, with plastic label

| | |
|-------------------|---------------|
| Width of box | 117 mm |
| Height of box | 82 mm |
| Depth of case | 440 mm |
| Width of case | 400 mm |
| Height of case | 129 mm |
| Case volume | 0.0227 cbm |
| Weight of product | 0.0074 kg |
| Weight of case | 4.52 kg |
| EAN of inner box | 4038917087315 |
| EAN case | 4038917005777 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

