



Product Specification

**S-Monovette® EDTA K3E, 2.6 ml, cap red, (LxØ):
65 x 13 mm, with plastic label pre-barcoded,
Pre-barcode with 8-digit unique number range
and 3-digit prefix**

Page 1



Product description

Order number	04.1901.088
Product description	S-Monovette® EDTA K3E, preparation: K3 EDTA, 2.6 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 011, label/print: white/red, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® EDTA K3E, 2.6 ml, cap red, (LxØ):
65 x 13 mm, with plastic label pre-barcoded,
Pre-barcode with 8-digit unique number range
and 3-digit prefix**

Page 2

Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	2.6 ml
Diameter	13 mm
Length including cap	81 mm
Length excluding cap	65 mm
Length including cap and plunger	102 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	42000
Depth of box	171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® EDTA K3E, 2.6 ml, cap red, (LxØ):
65 x 13 mm, with plastic label pre-barcoded,
Pre-barcode with 8-digit unique number range
and 3-digit prefix**

Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.004 kg
Weight of case	2.69 kg
EAN of inner box	4038917947947
EAN case	4038917947961

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

