



Product Specification

S-Monovette® Fluoride/EDTA FE, 2.6 ml, cap yellow, (LxØ): 65 x 13 mm, with plastic label

Page 1



Product description

Order number	04.1903
Product description	S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 2.6 ml, membrane screw cap, cap yellow, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label, label/print: transparent/yellow, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Fluoride + EDTA
Preparation concentration	1.0 mg/ml blood + 1.2 mg/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/yellow
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Fluoride/EDTA FE, 2.6 ml, cap yellow, (LxØ): 65 x 13 mm, with plastic label

Page 2

Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

Size

Sample volume	2.6 ml
Diameter	13 mm
Length including cap	81 mm
Length excluding cap	65 mm
Length including cap and plunger	102 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	42000
Depth of box	171 mm
Width of box	105 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Fluoride/EDTA FE, 2.6 ml, cap yellow, (LxØ): 65 x 13 mm, with plastic label

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.004 kg
Weight of case	2.69 kg
EAN of inner box	4038917087971
EAN case	4038917006132

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

