



Product Specification
S-Monovette® Serum CAT, 2.6 ml, cap white,
(LxØ): 65 x 13 mm, with plastic label

Page 1



Product description

Order number	04.1904
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 2.6 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label, label/print: transparent/white, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/white
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 2.6 ml, cap white,
(LxØ): 65 x 13 mm, with plastic label

Page 3

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0043 kg
Weight of case	2.85 kg
EAN of inner box	4038917087988
EAN case	4038917006149

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT