



Product Specification
S-Monovette® Lithium heparin LH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label

Page 1



Product description

Order number	04.1906.100
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 2.6 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 65 x 13 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	20 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label

Page 2

Liquid-tight Connection yes
membrane screw cap

Size

Sample volume 2.6 ml
Diameter 13 mm
Length including cap 81 mm
Length excluding cap 65 mm
Length including cap and plunger 102 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging case
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 42000
Depth of box 171 mm
Width of box 105 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0043 kg
Weight of case	2.86 kg
EAN of inner box	4038917134279
EAN case	4038917025584

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

