



Product Specification

S-Monovette® Homocysteine HCY Gel, 2.7 ml, cap light grey, (LxØ): 75 x 13 mm, with paper label

Page 1



Product description

Order number	04.1908.001
Product description	S-Monovette® Homocysteine HCY Gel, preparation: Homocysteine stabiliser, 2.7 ml, membrane screw cap, cap light grey, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: brown/grey, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Homocysteine stabiliser
Label/ Print	with paper label
Colour of print/label	brown/grey
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® Homocysteine HCY Gel, 2.7 ml, cap light grey, (LxØ): 75 x 13 mm, with paper label

Page 2

Size

Sample volume	2.7 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	light grey
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	171 mm
Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**S-Monovette® Homocysteine HCY Gel, 2.7 ml, cap
light grey, (LxØ): 75 x 13 mm, with paper label**

Page 3

Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0052 kg
Weight of case	3.36 kg
EAN of inner box	4038917163064
EAN case	4038917163057

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT