



Product Specification

S-Monovette® PFA, Citrate 9NC 0.129 mol/l 3.8% buffered, 3.8 ml, cap light blue, (LxØ): 65 x 13 mm, with paper label pre-barcoded

Page 1

Product description

Order number	04.1910.088
Product description	S-Monovette® PFA, Citrate 9NC 0.129 mol/l 3.8% buffered, preparation: Trisodium citrate/citric acid buffer solution, 3.8 ml, membrane screw cap, cap light blue, colour code ISO, (LxØ) without cap: 65 x 13 mm, with paper label pre-barcoded, label/print: light blue, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Trisodium citrate/citric acid buffer solution
Preparation concentration	0.129 mol/l pH 5.5
Colour code	ISO
Label/ Print	with paper label pre-barcoded
Colour of print/label	light blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	3.8 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	81 mm
Length excluding cap	65 mm
Length including cap and plunger	102 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	light blue
Membrane material	rubber
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® PFA, Citrate 9NC 0.129 mol/l 3.8% buffered, 3.8 ml, cap light blue, (LxØ): 65 x 13 mm, with paper label pre-barcoded

Page 2

Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	42000
Depth of box	171 mm
Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0045 kg
Weight of case	2.97 kg
EAN of inner box	4066916034092
EAN case	4066916034108

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT